

New Product Offering  
**Lucitone® FRS™ Flexible Partial Denture**  
By **Dentsply**

Dear Clinician,

We strive to keep you abreast of the latest advancements in dental technology by providing the best treatment plans and the most innovative products available.

You can now prescribe what we consider the finest Removable Partial Denture we have experienced in our combined 20-year history. Lucitone® FRS™ – Flexible Dental Resin offers all that your patients need in removable partial dentures:


- | <b>Aesthetics</b>  | <b>Comfort</b>  | <b>Durability</b>   |
|--|---|---|
| <ul style="list-style-type: none"><li>• Invisible Clear Clasps</li><li>• Four Lucitone 199® Shades and Clear</li></ul> | <ul style="list-style-type: none"><li>• True F-L-E-X-ability</li><li>• Monomer and Metal-Free</li></ul> | <ul style="list-style-type: none"><li>• Higher tensile &amp; flexural strength than other flexible acrylics</li></ul> |

**FRS is your most flexible option**

Lucitone FRS is a temporary thermoplastic injectable resin that offers superior flexibility, transparency for optimal aesthetics and stress crack resistance. Lucitone FRS material delivers first-in-class removable partial dentures. The partials are very lightweight and comfortable for your patients. Since the clasps on the partials are clear or tissue tone, they seem to disappear when in the mouth, making this an excellent alternative to a cast partial with clasps. And FRS delivers proven accuracy the first time, eliminating costly remakes.

Prescribe Lucitone FRS Flexible Partial Denture this week and see for yourself what a difference it makes. We ensure your satisfaction with a risk-free, money back offer .

Thank you for your time and attention,



Steve Leonetti  
President

"

"

....."; 58'UG'Cpngp{"Utggv."Uwkg'3."Rqtvrpf."QT"; 9463/3522".....": 9904: 507573

"

.....y y y @uwtgf f gpcmd@qo

"

"